



## PORTABLE OXYGEN & BACK-UP TANKS ON CONTINUOUS FLOW

Tank Size & Total Liters	<u>2 LPM</u>	<u>3 LPM</u>	<u>4 LPM</u>	<u>5 LPM</u>	<u>6 LPM</u>	<u>7 LPM</u>	<u>8 LPM</u>	<u>9 LPM</u>	<u>10 LPM</u>
B- Tank (M-6) 167 Liters	<b>1 hour 24 mins</b>	<b>55 mins</b>	<b>42 mins</b>	<b>33 mins</b>	<b>27 mins</b>	<b>23 mins</b>	<b>20 mins</b>	<b>18 mins</b>	<b>16 mins</b>
C- Tank 250 Liters	<b>2 hours 5 mins</b>	<b>1 hour 23 mins</b>	<b>1 hour 2 mins</b>	<b>50 mins</b>	<b>41 mins</b>	<b>35 mins</b>	<b>31 mins</b>	<b>27 mins</b>	<b>25 mins</b>
D-Tank 425 Liters	<b>3 hours 32 mins</b>	<b>2 hours 21 mins</b>	<b>1 hour 46 mins</b>	<b>1 hour 25 mins</b>	<b>1 hour 10 mins</b>	<b>1 hour</b>	<b>53 mins</b>	<b>47 mins</b>	<b>42 mins</b>
E-Tank 679 Liters	<b>5 hours 40 mins</b>	<b>3 hours 46 mins</b>	<b>2 hours 50 mins</b>	<b>2 hours 15 mins</b>	<b>1 hour 53 mins</b>	<b>1 hour 37 mins</b>	<b>1 hour 24 mins</b>	<b>1 hour 15 mins</b>	<b>1 hour 7 mins</b>

\*Note: If conserving device is used, supply is multiplied by a factor of 4-5

**445 Veit Rd, Suite C. Huntingdon Valley, PA 19006 PHONE: (215) 992-7770 FAX: (215) 992-7782**